

# QUICK GUIDE: CHARGER



1. Charging Coil
2. Charger Battery Level Indicator
3. Antenna Strength Indicator
4. Stimulation OFF Button
5. Charge Start Button
6. Power Adapter Port
7. IPG Battery Level Indicator



## PREPARING FOR CHARGING

1. Place the center of the circular Charging Coil over the top of your implanted pulse generator (IPG).

- Figure 1 shows that the top of the IPG (in yellow box) is at the center of the Charging Coil.
- Figure 2 shows an IPG implanted just below the belt-line. Your IPG may be in a different location, such as on your front side, in the abdomen (not shown).




Figure 1





Figure 2

## INITIATING CHARGING

2. Press  and hold until beeping starts, then release.

- Beeping and blinking of green lights will continue until the Charger links with the IPG. If a link with the IPG is not made within 10-12 seconds, the Charger will sound one long beep and shut off automatically. Reposition the Charging Coil and press  again.
- Once lights show on the Antenna Strength Indicator, a link has been made. More lights mean a better link, making it less likely that charging will be interrupted by small movements.
- If the link is suboptimal or lost, the Antenna Strength Indicator will show only the small orange bar and beep intermittently.

3. If the link is not made, suboptimal, or lost, reposition the blue Charging Coil *slowly* and repeat Steps 1 & 2.

## CHARGING

4. Once the link is made, you will see lights cascading upwards on both the Charger Battery Level Indicator and the IPG Battery Level Indicator, indicating charging is in progress.

- If during this process the link is lost, the Charger will beep as it searches for the IPG link. Reposition the blue Charging Coil and repeat Steps 1 & 2.

## COMPLETING CHARGING

5. Once the charging is completed, one long beep will sound and the lights on the Charger will turn OFF. Check the IPG Battery Level to make sure the battery has been fully charged.

6. Use the Remote Control to check that the IPG is turned ON.

7. Plug the Charger into the wall socket to ensure that the Charger will be ready for the next charging session.


- To plug the Charger into the wall socket, connect the provided Power Adapter into the Charger's Power Adapter Port first. Then plug the Power Adapter into the wall socket.
- There is no harm in continuously charging the Charger. The battery will not be damaged.




## USING THE CHARGING BELT

Once you know where to place the Charging Coil to achieve a strong link with the IPG, use the Charger Belt to secure the Charging Coil's location on your body *before charging*.

1. Place the blue Charging Coil into the mesh pouch on the Charger Belt and securely attach belt around the waist.
2. Initiate charging and move the Charging Coil until it links with the IPG.
3. If you want to be more mobile during charging, insert the Charger into the holster. Then, clip the holster to your belt or trousers. The holster should not be placed on the Charger Belt as it could displace the Charging Coil.



## TIPS

- Do not press  unless you want to turn OFF the therapy. In order to turn the therapy ON, use your Remote Control.
- Charging times may vary for each session, depending on factors such as recharging frequency and stimulation settings.
- If the Charging Coil has a suboptimal link as shown by the Antenna Strength Indicator, reposition the Charging Coil *slowly* and wait for 3 seconds. The Antenna Strength Indicator display may take up to 3 seconds to be updated for the new Charging Coil position. Repeat as necessary.
- If the Charging Coil loses the link with the IPG, a long beep will sound and the Charger will turn OFF automatically.
- If you want to turn OFF the Charger or discontinue charging, move the Charging Coil away from the IPG. The Charger will turn OFF after 10 seconds.
- If charging time is shorter than expected, then your IPG may have been turned OFF. Check to make sure your IPG is turned ON with the Remote Control.

## IMPORTANT REMINDERS

- Do not charge the IPG while the Charger is plugged into the wall. Always remove the Power Adapter before pressing the Charge Start Button.
- Recharge the IPG regularly to develop an easy-to-remember routine.
- While charging, the Charging Coil may become warm. If pain or discomfort is felt, cease charging and contact your doctor or Nevro representative.
- After recharging the IPG, plug the Charger back into the wall so that the Charger will be ready for your next recharging session.
- Use the Remote Control after charging to verify that the therapy is ON.
- If charging times dramatically increase, please notify your Nevro representative.

**For full instructions on this device, warnings, and precautions, please refer to the Patient Manual.**

FOR MORE INFORMATION  
ON HF10 THERAPY,  
TALK TO YOUR NEVRO  
REPRESENTATIVE.

All patients do not respond the same way to spinal cord stimulation (SCS) and experiences may vary. Patients should consult a physician to understand the potential benefits and risks of treatment with SCS.

**CE**  
0086 Effective  
4 May 2010

**USA**

**EC REP**

MDSS GMBH  
Schiffgraben 41  
D-30175 Hannover,  
Germany

Australian Sponsor:  
Emergo Asia Pacific Pty Ltd  
201 Sussex Street, Darling Park,  
Tower II, Level 20  
Sydney, NSW 2000  
Australia



**Nevro Corp.**

1800 Bridge Parkway  
Redwood City, CA 94065 USA  
[www.nevro.com](http://www.nevro.com)  
Tel. 1.650.251.0005  
Toll-free 1.888.956.3876

Senza, HF10, Nevro, and the Nevro logo are trademarks of Nevro Corp.  
©2015 Nevro Corp. All Rights Reserved.  
12125 Rev. B